



**Gronamic®  
Sport**

# Gronamic Sport High N

## High nitrogen levels for sustained growth

16-3-6 +2MgO +2CaO +0.5Fe  
(16-1.3-5 +1.2Mg +1.4Ca +Fe)

Gronamic Sport High N delivers sustained nutrition to the turf plant through the growing season. The supply of all key macronutrients and iron allows for a healthy plant, great turf colour and resilience from wear.

The balance of mineral and organic components provides the benefits from both approaches of turf nutrition. Supplies carbon to the soil/rootzone.

## Application notes

- ✓ Apply to dry foliage
- ✓ Water in if no rain falls within 1-2 days after application
- ✓ Irrigation or rainfall will assist dispersion and minimize risk on close mown surfaces
- ✗ Avoid applications during frosty or drought conditions
- ✗ Delay cutting or grooming for 2 days after application to allow product dispersal
- ✗ Product contains iron so if spilt on pavement, concrete or decorative surfaces brush off immediately as there is risk of staining



## Benefits

- ✓ Ideal through main growth periods of the season
- ✓ High nitrogen levels for sustained growth
- ✓ Combination of ammonium, urea and organic nitrogen sources
- ✓ Polyhalite mineral included to supply CaO, MgO & K for healthy growth
- ✓ Organic component supplies carbon to the soil

## Product characteristics

Formulation type:	Standard granule
Technology:	Slow release
Granule size:	2mm–4mm
Pack size:	25kg
Pack coverage:	625-833m <sup>2</sup>
Granule dispersal:	Less than 4 days
Longevity:	6-8 weeks
Typical turf response:	3-7 days (depending on environmental factors)

## Areas of use

Tees	★
Fairways	★
Stadia Pitches	★
Sports Pitches	★
Outfields	★
Racecourses	★
Lawns	★

## Application information

Application rate: **30-40g/m<sup>2</sup>**

## Gronamic Sport High N

### Product name

oxide analysis  
(elemental analysis)

### Product notes

### Recommended period of use

**Gronamic Sport High N**  
16-3-6 +2MgO + 2CaO +0.5Fe  
(16-1.3-5 +1.2Mg +1.4Ca +Fe)

Balanced mineral and organic nutrients provide available and slow release nutrition.

J F M A M J J A S O N D

## Spreader settings

### SR-2000 and Accupro 2000 spreader

Rate	Cone setting	Spread width	Single pass	Double pass
30g/m <sup>2</sup>	6	4.0m	V ½	O
35g/m <sup>2</sup>	6	4.0m	n/a	P
40g/m <sup>2</sup>	6	4.0m	n/a	Q